

The background features a large, stylized circular graphic on the right side, composed of overlapping geometric shapes in shades of orange, red, and white. A central portion of this circle is filled with a blurred, horizontal light trail in warm tones. The overall design is clean and modern, with a focus on the brand name and tagline.

Garrison[®]

Dental Solutions

Putting ideas
into practice

Product Catalog **2022**

☘ table of contents

MATRIX SYSTEMS

sectional matrix systems	Composi-Tight® 3D Fusion™	4-5
	NiTin™	8-9
	Matrix Rings	10-11
retainerless Tofflemire-style	Reel Matrix™	30-31
	Pinch™ Matrix	34
traditional Tofflemire-style	Slick Bands™	32-33
anterior matrix systems	Fusion™	6-7
	VariStrip™	36
cervical matrix systems	Blue View™ Cervical Matrices	37

WEDGES

plastic	Composi-Tight® 3D Fusion™ Wedge	16
	A+Wedge®	17
	WedgeWands®	17
	G-Wedge™	17
wood	SoftWedge™	17
pre-wedge & shield	FenderWedge®	19

IMPORTANT ACCESSORIES

interproximal finishers	FitStrip™	20-23
polishing systems	Rally™	24-25
	eZr™	26-27
forceps	Composi-Tight® Matrix Forceps	38
	NiTin™ Forceps	38
	Composi-Tight® 3D Fusion™ Forceps	39
	Ring Placement Forceps	39
	Universal Placement Forceps	39
specialty	PerForm™ Contact Former	35
composite instruments	Multi-Function Non-Stick	40
	Standard Non-Stick	41


INFORMATION PAGES

Composi-Tight® 3D Fusion™	2-3
Band Comparison Chart	13-15
Cases	12, 28, 29, 34
CE Information	18
Matrix Selection Guide	42-43
Order Form	44
Ordering Information	45
Guarantee	back cover

Only a sectional matrix can deliver results like this!

Tofflemire-style System

- ❌ Fails to restore proximal anatomy
- ❌ Thin contact at the marginal ridge
- ❌ Large food trap below
- ❌ Increased likelihood of fracture, occlusal interference, recurrent caries and periodontal disease



Composit-Tight® 3D Fusion™ Sectional Matrix System

- ✓ Operator-friendly retaining system
- ✓ Naturally contoured bands
- ✓ Anatomically correct contacts
- ✓ Contacts at the height of contour
- ✓ Rings produce optimal tooth separation for excellent, tight contacts

Only Composit-Tight® 3D Fusion™ can make the technique this intuitive and frustration free!



1 Place the matrix band:

Dead-soft three dimensionally contoured matrices have placement tabs for easier occlusal placement and Garrison's revolutionary non-stick coating for hassle-free removal.



2 Wedge firmly:

Composit-Tight® 3D Fusion™ Ultra-Adaptive wedges glide in but won't back out. Their soft, flexible exterior readily adapts to interproximal irregularities for a superior gingival seal.



3 Place the ring:

All three Composit-Tight® 3D Fusion™ separator rings will fit right over the wedge for easier placement. With Ultra-Grip™ retention extensions, they will stay firmly in place from the distal of the canine to the widest of molar preparations.

Composi-Tight® 3D Fusion™ Cases

Great adaptation to the teeth eliminates the need for adjustment with polishers in the proximal areas.



Images courtesy of Dr. Clarence Tam, Auckland, NZ

Wide preps can be handled in one simple step when you have the right (FX600 green) ring.



Images courtesy of Dr. Nikhil Sethi, BDS (Lon) Msc (Lon), Birmingham, UK

Back-to-back requires a shift in technique—we can show you how.*



Images courtesy of Dr. Jeffrey Kelley, DDS, Zeeland, MI

Quadrants can be a challenge. Ease your concerns with tips and tricks from the pros.



Images courtesy of Dr. Mike Malone, DDS, Lafayette, LA

*To see full writeups, cases and reviews go to www.garrisondental.com/learning
Do you have a case you want to share? Email it to us at contact@garrisondental.com.

The Award Winning Composi-Tight® 3D Fusion™ Sectional Matrix System



TESTIMONIAL

"3D rings are the real topic of Garrison's systems. The "v" shape of the ring that fits in the interproximal area allows a good fitting between the cavity margins and the matrix in the buccal and palatal walls. This results in an easier positioning of the composite masses close to the cavity margins and also the final remodeling usually necessary at the time of removal of the matrix will be very minimal. The rings also permit a divergence of the interproximal dental elements, and this causes a great point of contact. Garrison systems make Class II restorations simpler and more predictable and also reduce the operating time of the finishes when the matrix is taken off."

Enrico Cogo, DDS
Legnago, Ferrara, and Goito, Italy

Composi-Tight® 3D Fusion™

More than 20 years of sectional matrix leadership and an incredible amount of feedback from users all over the globe were poured into the design. All this hard work has resulted in separator rings with significantly improved performance and ease of use.



Tight contacts - guaranteed!


Superior tooth separation force through advanced engineering and material selection guarantees perfect contacts every time.

Simplified technique!

Composi-Tight® 3D Fusion™ is packed with features like easy ring stacking, over-the-wedge placement and Ultra-Grip™ anti-ring slip extensions that make system use intuitive and frustration free.

Ideal anatomy!

Soft-Face™ silicone ring tips with enhanced marginal ridge anatomy hug the pre-contoured dead-soft matrices to the tooth for flash free, perfectly contoured restorations.



Wider indications, wider smiles!

Dentistry's first Wide Preparation separator ring (green) makes Composi-Tight® 3D Fusion™ the most versatile and user-friendly sectional matrix system available. The Wide Prep ring greatly simplifies what until now has been one of the most challenging posterior composite restorations. Now you can restore these big preps with confidence knowing that you will achieve excellent contact and contour without having to take extra steps.

Composi-Tight® 3D Fusion™

Sectional Matrix System

A system designed to handle even the toughest cases



Composi-Tight® 3D Fusion™ Bands # .0016" thick	FX100	FX150	FX175	FX200	FX300
Non-stick dead soft stainless steel	bicuspids 4.0 mm w/ext 4.4 mm	bicuspids with deep cervical preparations 4.1 mm w/ext 6.0 mm	bicuspids and small molars 5.0 mm w/ext 5.6 mm	molars 5.4 mm w/ext 6.6 mm	molars with deep cervical preparations 6.0 mm w/ext 8.7 mm

Kit

- 1 Short Ring (blue)
- 1 Tall Ring (orange)
- 1 Wide Prep Ring (green)
- 70 Assorted Full Curve Non-Stick Matrix Bands
- 80 Assorted Ultra-Adaptive Wedges
- 4 Assorted FenderWedges® (pg 19)
- 12 Assorted Rally™ Polishers (pg 24-25)
- 1 Ring Placement Forceps (pg 39)
- 1 BFA-1 Band Placement Forceps (pg 38)
- 1 TN009 Multi-Function Instrument (pg 40)

Item

Item	Price
FXR-KFF-11	\$779
FXR-KFF-10 (BFA-1 and TN009 not in kit)	\$569
FXR-KFF-14 (TN009 not in kit)	\$689

Composi-Tight® 3D Fusion™ Rings

- Short Ring (blue)
- Tall Ring (orange)
- Wide Prep Ring (green)
- 3-Ring Kit (1 of each ring)

Item	Qty.	Price
FX400	2	\$242
FX500	2	\$242
FX600	2	\$242
FXR01	3	\$354

Composi-Tight® 3D Fusion™ Bands

- Bicuspids
- Bicuspids w/ deep cervical preps
- Bicuspids and small molars
- Molars
- Molars w/ deep cervical preps

	Standard (-M)		Bulk	
Item	Qty.	Price	Qty.	Price
Bicuspids	50	\$83	100	\$146
Bicuspids w/ deep cervical preps	30	\$83	60	\$146
Bicuspids and small molars	50	\$83	100	\$146
Molars	50	\$83	100	\$146
Molars w/ deep cervical preps	30	\$83	60	\$146

See cases and videos online!





Case using the Fusion™ Anterior Matrix System

Dr. Clarence "Clarie" Tam (Auckland, New Zealand) demonstrates her technique for efficient Class III restorations while bonding to multiple substrates.



Figure 1.
The final preparation form involved a substrate of both feldspathic porcelain/leucite-reinforced glass ceramic as well as dentin and cementum.



Figure 2.
The distal aspect of #8 was isolated with a Fusion™ Anterior sectional matrix and a curved orange wedge secured to ensure linguoproximal adaptation and thus a minimal of flash or overhangs of restorative material.



Figures 3 and 4.
The mesial of #8 was restored next, followed by the MP of #7.




Figure 5.
Completed restorations.

Fusion™ Anterior Matrix System

The **Fusion™ Anterior Matrix System** from Garrison® gives you the materials and techniques to restore many types of anterior lesions with confidence.

For years, our customers have shown us how they use posterior metal matrices to restore anterior teeth. When placed vertically, what was the buccal-lingual curvature of the band creates ideal gingival-incisal curvature for an anterior tooth. We worked with these customers to design a metal matrix band that takes this technique to a whole new level of refinement.

Add to that the intelligent design of the new **Fusion™ Anterior Wedges** and you get a complete system. These radically curved wedges hug the bands to the tooth all the way from the facial aspect around to the lingual, helping to free up your hands. Composite placement and sculpting are simplified, allowing you to easily express your artistic nature restoring beautiful smiles for your patients.



\$188

ANK01 Contains:

- 50 short bands
- 50 tall bands
- 25 extra small wedges
- 25 small wedges
- 25 medium wedges
- 25 large wedges



Fusion™ Anterior Bands

Refills	Item	Standard		Bulk	
		Qty.	Price	Qty.	Price
Short	AN100-M	50	\$82	AN100	100 \$144
Tall	AN200-M	50	\$82	AN200	100 \$144

Fusion™ Anterior Wedges

Description	Item	Qty.	Price	Item	Bulk	
					Qty.	Price
Fusion Anterior Wedge Kit	ANK4	400	\$181			
Fusion Anterior Wedge Kit	ANK4-M	200	\$95			

Refills	Item	Standard		Bulk	
		Qty.	Price	Qty.	Price
X-Small Yellow	ANYL-M	50	\$28	ANYL	100 \$51
Small Blue	ANBL-M	50	\$28	ANBL	100 \$51
Medium Orange	ANOR-M	50	\$28	ANOR	100 \$51
Large Green	ANGR-M	50	\$28	ANGR	100 \$51



NiTinol Sectional Matrix System

NiTinol + PEEK for the ultimate nickel titanium ring performance

As the name implies, the NiTin™ sectional matrix system separator rings are made from a special formulation of NiTinol nickel titanium for increased strength and resiliency. Rather than typical material used by other brands, NiTin™ is made with drawn-wire Nitinol which produces superior spring characteristics. Add in the PEEK reinforcement on the bow of the ring and you've got a system that produces exceptional contacts for hundreds and hundreds of uses.


Outperforms all other higher priced competitors' rings!



Smaller, less "bulky" design fits well in tight spots and is still stackable for MODs.

Drawn-wire nickel titanium for the ultimate in long life!


Proprietary heat-set processing increases strength and resiliency.



PEEK (polyether ether ketone) super-polymer reinforcement reduces ring stretch.

Ultra-durable PEEK tips are the strongest, longest lasting polymer tips on any ring.

Precisely engineered tip geometry allows NiTin™ to grip where others pop-off all while reducing flash.





NiTinol Sectional Matrix System



NiTIn™ Bands # .0016" thick					
Non-stick dead soft stainless steel	NTG100 bicuspids 4.0 mm w/ext 4.4 mm	NTG150 bicuspids with deep cervical preparations 4.1 mm w/ext 6.0 mm	NTG175 bicuspids and small molars 5.0 mm w/ext 5.6 mm	NTG200 molars 5.4 mm w/ext 6.6 mm	NTG300 molars with deep cervical preparations 6.0 mm w/ext 8.7 mm

The Complete NiTin™ System



NTG-KMG-01

- Kit**
- 1 Long NiTin™ Ring
 - 1 Standard NiTin™ Ring
 - 100 Assorted NiTin™ Matrices
 - 100 Assorted G-Wedges™
 - 1 NiTin™ Ring Placement Forceps
 - 1 NiTin™ Band Forceps

Item **Price**
NTG-KMG-01 \$455

- NiTIn™ Rings**
- Long Ring (black)
 - Standard Ring (white)

Item **Qty.** **Price**
NTG400 2 \$211
NTG500 2 \$211

- NiTIn™ Bands**
- Bicuspids
 - Bicuspids w/ deep cervical preps
 - Bicuspids and small molars
 - Molars
 - Molars w/ deep cervical preps

Item **Qty.** **Price**
NTG100 50 \$66
NTG150 50 \$66
NTG175 50 \$66
NTG200 50 \$66
NTG300 50 \$66

Composi-Tight® 3DFusion™

Sectional Matrix System



Simply the best!



Hugs the curves – Soft-Face™ silicone tips with enhanced marginal ridge anatomy hug the matrix band to the tooth for flash free restorations



Stack it up – Tru-Stack™ technology allows placement of the orange ring (tall) over the blue ring (short) for easier MOD and multiple tooth restorations



It just keeps going...
Enhanced durability for extended ring life and greater tooth separation

Get a grip! All Composi-Tight® 3D Fusion™ separator rings incorporate Ultra-Grip™ retention extensions that anchor them securely to virtually any tooth. These rings stay where you place them!



Composi-Tight. 3DFusion™

Soft-Face™ silicone tips with enhanced marginal ridge anatomy

\$242/pk

		
FX400 Short Ring 2/pk	FX500 Tall Ring 2/pk	FX600 Wide Prep Ring 2/pk

Composi-Tight. 3D XR ❖❖

The tenacious grip, along with enhanced Soft-Face™ technology and Dynamic Tip Angle™

\$236/pk



3DXR
Soft-Face™ Extra Retention Rings
2/pk

Composi-Tight. 3D ❖❖

The Soft-Face™ 3D Ring provides maximum adaptation to the tooth

\$231/pk



3D500
Soft-Face™ Rings
2/pk

Composi-Tight. 3D Clear ❖❖

The clear Soft-Face™ allows the curing light to be applied from both the buccal and lingual surfaces

\$231/pk



3D600
Clear Soft-Face™ Rings
2/pk

NiTiⁿ

Nitinol + Peek for the ultimate nickel titanium performance

\$211/pk



NTG400
Nitinol™ Long
2/pk



NTG500
Nitinol™ Standard
2/pk

Composi-TightGold. ❖❖

2x stronger than our Original G-Rings™

\$101/pk



AU400
Standard Length Tines
3/pk



AU500
Long Length Tines
3/pk

Composi-Tight. ❖❖

The original G-Ring™ that started it all—made of polished stainless-steel

\$91/pk



G100
Standard Length Tines
3/pk



G200
Long Length Tines
3/pk



Case using Composi-Tight® 3D Fusion™ Firm Matrix Bands

Dr. Anand Narvekar (Mumbai, India) utilized Garrison® Composi-Tight® 3D Fusion™ Firm Bands for deep margin elevation in this challenging onlay case.



Figure 1. FXH300 large molar band with extension was placed. The wedge was wrapped with multiple layers of teflon tape to sufficiently increase its size to allow it to seal the cervical portion of the matrix band in this extremely wide interproximal space.



Figure 2. Cervical area was built up with multiple layers of composite.



Figure 3. The process was repeated for the adjacent molar. The stiff nature of Composi-Tight® 3D Fusion™ Firm Bands allows apical pressure to be applied to the band without deformation of the matrix. It can be seated deeper than dead soft matrices for these types of procedures.



Figure 4. One-week post-bonding of onlays.

Composi-Tight® 3D Fusion™ Firm Matrix Bands



Composi-Tight® 3D Fusion™ Firm

matrix bands with eZ-Place™ technology maintain their shape during placement making them *significantly easier* to insert!

eZ-Place™ technology combines specially processed stainless steel with a proprietary matrix band forming process to yield Fusion Firm's remarkable handling characteristics. They resist deformation nearly 70% better than typical dead-soft bands yet are pre-contoured and super thin at 0.0015".

Fusion Firm:

- Nearly 70% better at resisting deformation!
- Super-thin at only 0.0015"
- Five ideal sizes!






Perfect for:

- Tight interproximal spaces!
- More conservative preps!
- First-time users and veterans alike!



0.0015" thick

FXH Series

	Bicuspids	Bicuspids with deep cervical preps	Bicuspids and small molars	Molars	Molars with deep cervical preps
Composi-Tight® 3D Fusion™ Firm	Composi-Tight® 3D Fusion™ Firm matrix bands with eZ-Place™ technology for easier placement.				
thickness 0.0015"					
height	FXH100 4.3 mm	FXH150 4.3 mm w/ext 5.9 mm	FXH175 5.3 mm	FXH200 6.1 mm	FXH300 6.1 mm w/ext 7.7 mm


FX Series

	Bicuspids	Bicuspids with deep cervical preps	Bicuspids and small molars	Molars	Molars with deep cervical preps
Composi-Tight® 3D Fusion™	Composi-Tight® 3D Fusion™ Full Curve Bands have non-stick coating and an occlusal placement tab.				
thickness 0.0016"					
height	FX100 4.0 mm w/ext 4.4 mm	FX150 4.1 mm w/ext 6.0 mm	FX175 5.0 mm w/ext 5.6 mm	FX200 5.4 mm w/ext 6.6 mm	FX300 6.0 mm w/ext 8.7 mm

Description	Item	Qty.	Kit	Standard (-M)		Bulk	
				Qty.	Price	Qty.	Price
Assorted Band Kit	FXHB04	300	\$355				
Assorted Band Mini Kit	FXHB05	150	\$185				

Description	Item	Qty.	Kit	Standard (-M)		Bulk	
				Qty.	Price	Qty.	Price
Assorted Band Kit	FXB01	420	\$646				
Assorted Band Mini Kit	FXB02	210	\$331				

- FX Series
- FXH Series
- NTG Series
- SXR Series
- B Series
- M - SM - AU Series



Band Comparison Chart



NTG Series





NiTin™ NiTin™ Full contour matrix bands for ideal anatomy.

thickness 0.0016"	NA					
height		NTG100 4.0 mm w/ext 4.4 mm	NTG150 4.1 mm w/ext 6.0 mm	NTG175 5.0 mm w/ext 5.6 mm	NTG200 5.4 mm w/ext 6.6 mm	NTG300 6.0 mm w/ext 8.7 mm

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
NTG100	25	\$37	50	\$66
NTG150	25	\$37	50	\$66
NTG175	25	\$37	50	\$66
NTG200	25	\$37	50	\$66
NTG300	25	\$37	50	\$66

SXR Series

Slick Bands™ XR Slick Bands™ XR Matrices have non-stick coating and Grab-Tabs™ for easier placement.

thickness 0.0016"	NA					
height		SXR100 4.6 mm	SXR150 3.8 mm w/ext 5.4 mm	SXR175 5.4 mm	SXR200 6.4 mm	SXR300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
SXR100	50	\$75	100	\$131
SXR150	30	\$75	60	\$131
SXR175	50	\$75	100	\$131
SXR200	50	\$75	100	\$131
SXR300	30	\$75	60	\$131

SM Series

Slick Bands™ Slick Bands™ Matrices have non-stick coating.

thickness 0.0016"						
height	SM050 3.2 mm w/ext 4.2 mm	SM100 4.6 mm	SM150 3.8 mm w/ext 5.4 mm	SM175 5.5 mm	SM200 6.4 mm	SM300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
SM050	NA	NA	60	\$131
SM100	50	\$75	100	\$131
SM150	30	\$75	60	\$131
SM175	50	\$75	100	\$131
SM200	50	\$75	100	\$131
SM300	30	\$75	60	\$131

To order standard packs, add a -M after the item #.

AU Series

	Pediatric	Bicuspid	Bicuspid with deep cervical preps	Bicuspid and small molars	Molars	Molars with deep cervical preps
Composi-Tight Gold*	The Gold bands are the same as the Composi-Tight® bands.					
thickness 0.0015"				NA		
height	AU050 3.2 mm w/ext 4.2 mm	AU100 4.6 mm	AU150 3.8 mm w/ext 5.4 mm	NA	AU200 6.4 mm	AU300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
AU050	25	\$62	50	\$95
AU100	50	\$62	100	\$95
AU150	25	\$62	50	\$95
AU200	50	\$62	100	\$95
AU300	25	\$62	50	\$95

M Series

	Pediatric	Bicuspid	Bicuspid with deep cervical preps	Bicuspid and small molars	Molars	Molars with deep cervical preps
Composi-Tight®	Composi-Tight® bands are contoured in 3 dimensions to provide an accurate contact at the tooth's natural height of contour.					
thickness 0.0015"						
height	M050 3.2 mm w/ext 4.2 mm	M100 4.6 mm	M150 3.8 mm w/ext 5.4 mm	M175 5.5 mm	M200 6.4 mm	M300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
M050	25	\$62	50	\$95
M100	50	\$62	100	\$95
M150	25	\$62	50	\$95
M175	50	\$62	100	\$95
M200	50	\$62	100	\$95
M300	25	\$62	50	\$95

B Series

	Pediatric	Bicuspid	Bicuspid with deep cervical preps	Bicuspid and small molars	Molars	Molars with deep cervical preps
Composi-Tight® Original	The original bands have a slightly shorter length.					
thickness 0.0015"				NA		
height	B050 3.2 mm w/ext 4.2 mm	B100 4.6 mm	B150 3.8 mm w/ext 5.4 mm	NA	B200 6.4 mm	B300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
B050	25	\$54	50	\$81
B100	50	\$54	100	\$81
B150	25	\$54	50	\$81
B200	50	\$54	100	\$81
B300	25	\$54	50	\$81

MC Series

	Pediatric	Bicuspid	Bicuspid with deep cervical preps	Bicuspid and small molars	Molars	Molars with deep cervical preps
Composi-Tight® 3D Clear	Clear bands are transparent and translucent to allow for trans-enamel polymerization without interference from metal bands.					
thickness 0.0020"				NA	NA	NA
height	NA	MC100 5.0 mm	NA	MC175 5.9 mm	MC200 6.8 mm	NA

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
MC100	50	\$75	100	\$131
MC175	50	\$75	100	\$131
MC200	50	\$75	100	\$131

To order standard packs, add a -M after the item #.

Composi-Tight® 3D Fusion™

Ultra Adaptive Wedges

Combining Soft-Face™ adaptive materials, a firm inner core and advanced mechanical features to produce a wedge that truly works. Every time.



Fins with serious function
Soft retentive fins smoothly fold down during wedge insertion and then spring back when clear of the interproximal space. Wedge back-out is a thing of the past.



A real softie Composi-Tight® 3D Fusion's Soft-Face™ over-mold allows the wedge to do what no other wedge can truly do - actually adapt to interproximal irregularities.

Tough on the inside

The firm inner skeleton allows for easy insertion and tooth separation just like a traditional wedge.

Composi-Tight® 3D Fusion™ Wedge

Description

Composi-Tight® 3D Fusion™ Wedge Kit
Composi-Tight® 3D Fusion™ Wedge Kit

Item	Qty.	Kit
FXK4	400	\$213
FXK4-M	200	\$112

Refills

X-Small Yellow	FXYL-M	50	\$33
Small Blue	FXBL-M	50	\$33
Medium Orange	FXOR-M	50	\$33
Large Green	FXGR-M	50	\$33


Item	Standard		Item	Bulk	
	Qty.	Price		Qty.	Price
FXYL	100	\$61	FXYL	100	\$61
FXBL	100	\$61	FXBL	100	\$61
FXOR	100	\$61	FXOR	100	\$61
FXGR	100	\$61	FXGR	100	\$61



A+Wedge®

Controls Gingival Bleeding

Controls gingival bleeding with an aluminum sulfate coating that has an astringent effect on the gingiva, improves clarity of restoration area and eliminates need to apply an astringent separately



Description
A+Wedge® Kit

Item	Qty.	Kit
GWAK4	200	\$117


Refills
X-Small Yellow
Small Blue
Medium Orange
Large Green

Item	Bulk	
	Qty.	Price
GWAYL	100	\$61
GWABL	100	\$61
GWAOR	100	\$61
GWAGR	100	\$61

G-Wedge™

Perfect wedge. Perfect contact.

If you prefer traditional wedge placement with either a hemostat or cotton pliers, G-Wedges give you the same great shape as WedgeWands®, just without the applicator wand.



Description
G-Wedge™ Kit

Item	Qty.	Kit
GWK4	400	\$158

Refills
X-Small Yellow
Small Blue
Medium Orange
Large Green

Item	Standard		Item	Bulk	
	Qty.	Price		Qty.	Price
GSWYL	100	\$46	GWYL	300	\$125
GSWBL	100	\$46	GWBL	300	\$125
GSWOR	100	\$46	GWOR	300	\$125
GSWGR	100	\$46	GWGR	300	\$125

WedgeWands®

WedgeWands® are designed to twist off the wand, but many prefer to leave it attached to help retract the cheek and for easy removal after the procedure is complete.



Description
WedgeWands® Kit

Item	Qty.	Kit
WK4	400	\$183

Refills
X-Small Yellow
Small Blue
Medium Orange
Large Green

Item	Standard		Item	Bulk	
	Qty.	Price		Qty.	Price
SWYL	100	\$55	WYL	300	\$147
SWBL	100	\$55	WBL	300	\$147
SWOR	100	\$55	WOR	300	\$147
SWGR	100	\$55	WGR	300	\$147

SoftWedge™

Superior adaptation and seal

Soft Wedges™ compress during insertion, then swell to help fill voids to ensure that your preparations stay clean and dry—an absolute must for today's bonded restorations!



Description
Soft Wedge™ Kit

Item	Qty.	Kit
WDK4	400	\$98

Refills
X-Small Yellow
Small Blue
Medium Orange
Large Green

Item	Standard		Item	Bulk	
	Qty.	Price		Qty.	Price
WSDYL	100	\$30	WDYL	300	\$81
WSDBL	100	\$30	WDBL	300	\$81
WSDOR	100	\$30	WDOR	300	\$81
WSDGR	100	\$30	WDGR	300	\$81

Great technique tips and product help are only a click away!

Join thousands of subscribers on our YouTube channel and get notified instantly when new content is posted.

The Garrison® YouTube channel is NOT a collection of infomercials. It's all about helpful information on a wide variety of topics. From short two-minute tips on how to reduce flash on Class IIs, to more in depth live patient case presentations, the Garrison® channel delivers great content.

Go to youtube.com/garrisdental or search Garrison Dental on youtube and subscribe today!



Speaking of great content...

Get CE credits direct from Garrison®!

Go to the Learning Center at garrisdental.com to see a list of pre-recorded and upcoming live webinars – many with available CE.

Register for an online account on our website to sign-up for CE credits and manage your CE courses.


FENDER WEDGE®

Interproximal Tooth Shields

Pre-separate and Protect

FenderWedge® prevents damage to the adjacent tooth and gingiva that can occur during Class II preparations.

FenderWedge® is a combination of a steel plate and a plastic wedge. Inserted into the interproximal space, it protects the tissue and separates the teeth, thus simplifying the application of a matrix.



FenderWedge®

Description
FenderWedge® Kit


Item **Qty.** **Kit**
TSK01 80 \$128

Refills	Item	Standard		Bulk	
		Qty.	Price	Qty.	Price
X-Small Purple	TSPR-M	30	\$56	TSPR	60 \$102
Small Orange	TSOR-M	30	\$56	TSOR	60 \$102
Medium Green	TSGR-M	30	\$56	TSGR	60 \$102
Large Yellow	TSYL-M	30	\$56	TSYL	60 \$102



FenderWedge is NOT a sectional matrix system nor does it replace the use of a matrix band.

Clinical Tip



Place FenderWedge® prior to preparation to pre-separate and protect the adjacent teeth.

The plastic wedge—not just the metal shield—should firmly engage the interproximal surfaces of the teeth. This will provide the optimum pre-separation force and hold the FenderWedge® securely in place. (Figures 1 and 2)

Figure 3 clearly shows marks from the bur along the shield. Much better to have those on the shield than on the adjacent tooth!

Once the preparation is complete, remove the FenderWedge® and place your matrix.

Dentistry and photos courtesy of Dr. Rob Ritter, Jupiter, Florida

Visit garrisdental.com/ learning to see Dr. Ritter's full case presentation.





Case using FitStrip™ Subgingival Trimmers and Finishers

Doctors R. Becciani and V. Attanasio of Italy present a host of cases highlighting the unique attributes of FitStrip™.



Figure 1. Interproximal reduction for orthodontic reasons. The color coded, double sided abrasive strips work in concert with the calibrated gauges to provide the appropriate interproximal space quickly and accurately for aligner systems.



Figure 3. Separate inadvertently bonded together adjacent teeth. The very thin serrated FitStrip™ (gray barrel) is incredibly useful for separating teeth without damaging the interproximal surfaces or losing the contact point.



Figure 2. Removal of excess bond and cement on indirect restorations. Excess bond or cement in the particularly delicate cervical area can be achieved without iatrogenic trauma to the patient.



Figure 4. Prepare the interproximal space for minimally invasive ceramic veneers. The FitStrip™ system allows for the calibrated and minimal opening of the interproximal space.

To see more of how doctors Becciani and Attanasio use FitStrip™ as well as FitStrip™ videos, go to our website at garrisondental.com and click on FitStrip™.




Flexible
Interproximal finishers with a
Twist

FitStrips are a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts, **Garrison® Dental Solutions**.

Do the Twist! Twist the color-coded barrel to adjust FitStrips' curvature and it's automatically locked in place.

One handy handle! FitStrip's simple, attachable handle makes a world of difference for both clinician ease-of-use and patient comfort.



Ready for a marathon! FitStrips are very durable and long lasting. They stay sharp and ready to go even after multiple autoclave cycles.

Starts bold – finishes clean

1. Medium grit for fast removal
2. Fine grit smooths things out
3. Super fine strips put the perfect finishing touch on things



Flexible
Interproximal finishers with a
Twist



The FitStrip™ kit that started it all!

Most popular!

- Primary Function: Contouring Composites and IPR up to .30mm.
- Durable diamond abrasive on just one side to protect adjacent teeth
- The FitStrip™ serrated saw is the PERFECT instrument to separate accidentally bonded teeth and reduce over-tight contacts

FitStrip™ Starter Kit

FPSK01 contains:, 2 handles, 2 serrated, 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)

\$184

4 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm)



Your Complete Solution for Interproximal Reduction (IPR)!

- Primary Function: IPR when .40mm and .50mm are needed, yet can still be used to contour composites
- IPR gauges for exact reduction
- Attachable handle for greater ease of use and patient comfort

Universal Kit

FPSK04 contains:, 2 handles, 1 IPR gauge, 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)

\$203

6 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm, .40 mm, .50 mm)



Subgingival FitStrips are the ideal solution for the cervical area!

- Primary Function: Use with Composites, both for contouring and accessing cervical areas to remove overhangs or excess cements.
- Only 2.5mm height keeps abrasive below the contact point
- Safety zone protects contacts during insertion
- Extra thin at .05mm

FitStrip™ Subgingival Finishing and Single Sided Kit

FPSK08 contains:, 2 handles, 10 Strips: 2 subgingival super fine, 2 subgingival fine, 1 subgingival medium, 1 serrated saw, 1 each single sided .08 mm, .10 mm, .13 mm, .18 mm

\$184



FitStrips – a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts, **Garrison Dental Solutions.**

Kits All kits come in a package of 10 FitStrips with 2 handles.

Description/Contents	Item	Price	Primary Use
Starter Kit for finishing/polishing 2 serrated 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm) 4 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm)	FPSK01 (shown on page 22)	\$184	Finishing composites
Universal Kit for IPR and finishing/polishing 1 IPR gauge 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm) 6 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm, .40 mm, .50 mm)	FPSK04 (shown on page 22)	\$203	Finishing composites and IPR
IPR Kit 1 IPR gauge 10 double sided (2 of each .11 mm, .21 mm, .30 mm, .40 mm, .50 mm)	FPSK05	\$203	IPR
Subgingival Finishing Kit 3 subgingival super fine, 3 subgingival fine, 3 subgingival medium, 1 serrated saw	FPSK07	\$184	Subgingival
Subgingival Finishing and Single Sided Kit 2 subgingival super fine, 2 subgingival fine, 1 subgingival medium, 1 serrated saw, 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)	FPSK08 (shown on page 22)	\$184	Subgingival and finishing composites

Refills	Single Sided	Double Sided	10 pk w/2 handles	4 ct packs (add - 4 to model #)
Super fine 15 Micron 	FPSXFSS 0.08 mm thick	FPSXFDS 0.11 mm thick	Either Refill \$184	\$88
Fine 30 Micron 	FPSFSS 0.10 mm thick	FPSFDS 0.15 mm thick	Either Refill \$184	\$88
Medium 46 Micron 	FPSMSS 0.13 mm thick	FPSMDS 0.21 mm thick	Either Refill \$184	\$88
Coarse 90 Micron 	FPSCSS 0.18 mm thick	FPSCDS 0.30 mm thick	Either Refill \$184	\$88
Serrated 	FPSSAW 0.05 mm thick		Refill \$184	\$88
Medium 46 Micron 		FPSTMDS 0.40 mm thick	Refill \$184	\$88
Medium 46 Micron 		FPSXTMDS 0.50 mm thick	Refill \$184	\$88
Super fine 15 Micron 	FPSGXFSS 0.08 mm thick	N/A	Refill \$184	\$88
Fine 30 Micron 	FPSGFSS 0.10 mm thick	N/A	Refill \$184	\$88
Medium 46 Micron 	FPSGMSS 0.13 mm thick	N/A	Refill \$184	\$88
FitStrip™ replacement handles (2-pack) FPSHANDLE-2			\$31	
IPR Gauge FPSGAUGE			\$41	



Rally™

Mini Polishers for Composite

Rally gets you to the finish line faster!

The composite finish line that is. One lightning fast step will put a beautiful satin finish on posterior restorations or add the gray high shine polisher for a lustrous gloss on facials and anteriors. And if you want to be quick off the starting line, use the blue coarse grit for rapid reduction and shaping. Rally Mini Polishers get you across the finish line and on to your next patient in record time!

Get to the point - Every time!

Single patient Rally means that you get a clean, fresh, properly shaped polisher every time.


A clean get-a-way!

You don't want to leave a trail of crumbs or spattered polishing paste on the track. Rally's diamond impregnated high shine final polishers won't crumble during use and don't require paste to create a beautiful, durable shine.

Mini and mighty!


Rally's small size and tapered shaft create a truly agile polisher capable of reaching the smallest anatomical features without blocking your visibility.






**Remove -
Choose the Blue
(silicon carbide)**

Polishers for fast removal and contouring



**Smooth -
Choose the Red
(aluminum oxide)**

Polishers for smoothing and silky shine



**Shine -
Choose the Gray
(diamond impregnated)**

Polishers for ultra high shine



Slow speed right angle latch shank. Use with water.

Rally™

Item	Description	Qty.	Price
FPCCK01	Composite Polisher Complete Kit	60 (10 each)	\$156

Item	Refills	Qty.	Price
FPCCK010	Blue coarse grit point composite polisher	30	\$61
FPCCK040	Blue coarse grit cup composite polisher	30	\$61
FPCCK010	Red medium grit point composite polisher	30	\$61
FPCCK040	Red medium grit cup composite polisher	30	\$61
FPCCK010	Gray fine grit point composite polisher	15	\$61
FPCCK040	Gray fine grit cup composite polisher	15	\$61



Intra-oral Adjustment/Polishing
for Zirconia and Lithium Disilicate

Adjust and Polish Zirconia and
Lithium Disilicate the **eZr** Way!

eZr Intra-oral Adjustment:



- **No chatter, chipping or microfracture** - Unlike diamonds bonded directly to a hard metal shaft, eZr's rubber substrate cushions the diamonds to prevent damage to the restoration.
- **No heat build-up** - The carefully engineered substrate does not transfer heat during use.
- **FAST!** - eZr's coarse grit blue grinders make your adjustments fast and convenient.

eZr Intra-oral Polishing:



- **Unparalleled smooth polish** - Don't take our word for it! - "The eZr system is the only one I've tried that polishes zirconia to smooth!" Todd Brower, DDS, Holland, Michigan.
- **Increased physical stability** - Highly polished zirconia and lithium disilicate is stable and non-abrasive to opposing dentition.
- **No paste required** - Forget the mess and expense of polishing pastes. eZr™ is all you need.
- **No water required** - No need to use water when polishing. Water is only required for high speed.

eZr™ is a three-step diamond polishing system for the processing of zirconium oxide and lithium disilicate. The grinders and polishers are ideal for adjusting and correcting chairside restorations such as inlays, onlays, veneers and crowns.



*Slow speed right
angle latch shank
Steam autoclaveable*




Start by gross adjusting and contouring using the blue coarse grit grinders. Use the flat end taper for occlusal/labial/lingual/buccal surfaces and the disc for interproximal surfaces.



Use the green medium grit polishers for smoothing and pre-polishing. Use the point for the occlusal grooves, the cup for labial/lingual/buccal, and the disc for interproximal surfaces.



Use the orange fine grit polishers for high gloss polishing. Use the point for the occlusal grooves, the cup for labial/lingual/buccal, and the disc for interproximal surfaces.

Item	Description	Qty.	Price
FPZK01	Complete system kit Gross removal, smoothing, pre-polishing, and high gloss polishing	2 grinders 6 polishers	\$202
FPZC030	 Blue coarse grit disc Gross removal	5/pk	\$162
FPZC010	 Blue coarse grit flat end taper Gross removal	5/pk	\$122
FPZM010 (not in kit)	 Green medium grit small point polisher Smoothing and pre- polishing	5/pk	\$49
FPZM020	 Green medium grit large point polisher Smoothing and pre- polishing	5/pk	\$76
FPZM030	 Green medium grit disc polisher Smoothing and pre- polishing	5/pk	\$76
FPZM040	 Green medium grit cup polisher Smoothing and pre- polishing	5/pk	\$76
FPZF010 (not in kit)	 Orange fine grit small point polisher High gloss polishing	5/pk	\$49
FPZF020	 Orange fine grit large point polisher High gloss polishing	5/pk	\$76
FPZF030	 Orange fine grit disc polisher High gloss polishing	5/pk	\$76
FPZF040	 Orange fine grit cup polisher High gloss polishing	5/pk	\$76

Cases using other matrix systems



Images courtesy of Dr. Julio Maya



Images courtesy of Dr. Dory Stutman

Reel Matrix™

TESTIMONIALS

"Using a matrix system like Reel Matrix can enable them (doctors) to quickly and precisely place anatomically correct Class II posterior composite restorations that demonstrate exceptional contacts."

Dr. Dory Stutman

"We produced the margin elevation band to facilitate the deep margin elevation technique which is becoming increasingly popular with adhesive dentistry. Properly performed, this technique allows for a conservative indirect restoration to be completed without additional invasive treatments. The use of this band and unique matrix shape allows the band to form a tight seal in deeper cervical areas where other bands would usually be unable to reach."

*Dr. Matthew A. Nejad, DDS
Beverly Hills, CA*



Images courtesy of Dr. Matthew A. Nejad, DDS, Beverly Hills, CA

Margin Elevation Band




Images courtesy of Dr. Jack D. Griffin, Jr., DMD, MAGD, Weldon Spring, MO

PerForm™

"To ensure a tight contact, a bulk fill flowable was placed in one proximal box at a time and the mylar held against the adjacent tooth with a PerForm™ proximal contact instrument. This ensures that the contact is tight and well-formed while light curing though it."

*Dr. Jack D. Griffin, Jr., DMD, MAGD
Weldon Spring, MO*

Cases using anteriors/cervicals



Images courtesy of Dr. Zunzarren, DDS



Images courtesy of Dr. Abdi Sameni, DDS

VariStrip™

How can you produce a natural curvature for anterior composites and eliminate black triangles? Check out Garrison's curved "Varistrips" to recreate the natural interproximal anatomy for anterior teeth.

- Helps avoid incorrect flat embrasures, black triangles and food traps that can be caused by old-fashioned straight mylars.
- No need for wedges on subgingival fillings, as the contours naturally close off the apical aspects for you.
- Excellent for class IVs, diastema closures, etc.
- Variable height allows you to slide each strip in until it precisely matches the tooth height before cutting off excess.
- The blue tint does not affect cure, but ends the annoyance of losing clear strips on patient bibs or bracket tables.

Blue View VariStrips—a must-have super-simple product for every dental office.

Dr. Michael Curtis, Bridgeport, CT, USA



Clinical Tip

After using your preferred bonding procedure, place the restorative material into the preparation in advance of using the cervical matrix.

Place the cervical edge of the matrix against the tooth (Fig. 1).

While pressing the gingival edge of the matrix firmly against the tooth, roll the matrix onto the tooth (Fig. 2). Excess material will be expressed from the occlusal edge of the matrix rather than subgingivally, greatly simplifying clean up.

Cure directly through the matrix.

ReelMatrix

Real Easy. Real Tight. Real Fast!

For all posterior restorations!

1 Real Easy:
The reels pop into the handle without any fuss! Press the reel into the handle and retract the grippers to lock the reel in place. Plus, if you use the regular matrices (yellow hubs) they can be easily pushed through contacts.

2 Real Tight:
Tight contacts and excellent contour are what Reel Matrix is all about. These thin matrices are pre-contoured in three dimensions and when combined with anatomical wedges and separator rings, the result is amazingly tight contacts at the natural height of contour.

3 Real Fast:
Pre-loaded reels mean you can just lock them in and go. The retainerless design greatly speeds up quadrant work while improving visibility and patient comfort. Add an FX500 separator ring and you virtually eliminate buccal and lingual flash.



See page 34
for Dr. Rosenfeld's
Margin Elevation case.

Margin Elevation Technique restores deep caries without the trauma.

Margin elevation requires a matrix that will properly adapt and isolate the subgingival area. This is not possible with traditional matrix materials (Figure A). The narrow Margin Elevation Band, with its enhanced conical design, is able to slide down subgingivally and properly seal this area (Figure B).




Figure A




Figure B

For the complete margin elevation technique, go to garrissdental.com



RM500H
Margin elevation



RM100D
5mm deadsoft



RM200D
6.3mm deadsoft



RM100H
5mm regular



RM200H
6.3mm regular



RM100B
5mm blue



RM200B
6.3mm blue

The Reel Matrix™ is the ideal delivery system for the Margin Elevation Band. Its extremely small size and light weight prevents the accidental dislodgement that can occur with long, cumbersome metal retainers.



Margin elevation band is .002" thin. Super-thin dead-soft bands are .0012" thin. Regular bands are .0015" thin and can be placed through existing contacts. A Blue View cure-through option at .0029" thin is available if you prefer trans-enamel curing.

ReelMatrix

Real Easy. Real Tight. Real Fast!

Item	Description	Price
RMK05	Kit Contains 80 matrices (20 RM500H, 10 ea. RM100D, RM200D, RM100H, RM200H, RM100B, RM200B) 1 tensioning instrument 1 PerForm™ instruments (set of 2)	\$485

Refills		Qty.	Price
RM500H	Margin elevation band	50	\$94
RM100D	5mm dead-soft matrices	50	\$85
RM200D	6.3mm dead-soft matrices	50	\$85
RM100H	5mm regular matrices	50	\$85
RM200H	6.3mm regular matrices	50	\$85
RM100B	5mm blue matrices	50	\$95
RM200B	6.3mm blue matrices	50	\$95
RMR100	5.9mm empty reel	100	\$55
RMR200	6.7mm empty reel	100	\$55
RMH-100	Tensioning instrument		\$247

IMPORTANT ACCESSORIES

PerForm™ Proximal Contact Instrument (included in RMK05 kit)	Page 35
---	---------

Great for core build-ups

SlickBands™

Tofflemire-style

Non-stick technology developed by Garrison®: Modern bonding agents can adhere to uncoated matrix bands- Slick Bands™ slide right out.

Reaches over to the adjacent tooth more easily to help you create broader, deeper contacts.

- Extremely wide preps are simplified.
- Improved visibility for core buildups.



The increased curvature of the all-new **Right-Curve matrix band** produces outstanding gingival seal and reaches over to adjacent teeth easier for broader, deeper contacts. This creates a significantly increased funnel-shaped cone when wrapped around a tooth. It's the right choice for many restorations.


Tighter gingival seal is beneficial for all restorations.

Traditional band with narrow 'funnel'.



Margin Elevation Band

The deep margin elevation technique is a useful non-surgical approach to treating deep localized margins.



A margin elevation matrix will give you the best gingival adaptation and contour for the deep margin areas. It is trimmed in height to allow improved adaptation to deep margin areas.

"We produced the margin elevation band to facilitate the deep margin elevation technique which is becoming increasingly popular with adhesive dentistry. Properly performed, this technique allows for a conservative indirect restoration to be completed without additional invasive treatments. The use of this band and unique matrix shape allows the band to form a tight seal in deeper cervical areas where other bands would usually be unable to reach."

*Matthew A. Nejad, DDS
Beverly Hills, CA*

Slick Bands™ Tofflemire-style Matrices



Clinical Tip: Use dead-soft green bands when contact is COMPLETELY broken with adjacent teeth for maximum burnishability. When contact is NOT broken with adjacent teeth, use regular gray bands.

IMPROVED Dead-soft Bands


with non-stick coating

Regular Bands


with non-stick coating




Right-Curve™ matrices provide a lower contact point and improved isolation. Excellent for all restorations.



#13 Pedodontic for short and primary teeth



#2 Sub-gingival for deep restorations



Tofflemire-style matrices are .0014" thick
Margin elevation matrices are .002" thick



See page 34
for Dr. Rosenfeld's
Margin Elevation case.



Slick Bands™

Description

Slick Bands Kit
300 Assorted Slick Bands
100 G-Wedges
PerForm instruments (1 set; shown on page 35)

Item

SMTK06

Kit

\$349

Dead Soft Bands

Item	Qty.	Price
#1 Universal Tofflemire	50	\$42
Right-Curve	50	\$54
#13 Pedodontic	50	\$42
#2 Sub-gingival	50	\$42

Standard (-M)		Bulk	
Qty.	Price	Qty.	Price
100	\$42	100	\$67
100	\$54	100	\$82
100	\$42	100	\$67
100	\$42	100	\$67

Regular Bands

Item	Qty.	Price
#1 Universal Tofflemire	50	\$38
Right-Curve	50	\$50
#13 Pedodontic	50	\$38
#2 Sub-gingival	50	\$38

Standard (-M)		Bulk	
Qty.	Price	Qty.	Price
100	\$38	100	\$60
100	\$50	100	\$79
100	\$38	100	\$60
100	\$38	100	\$60

Regular Bands

Item	Qty.	Price
Margin elevation band	50	\$59

Standard (-M)		Bulk	
Qty.	Price	Qty.	Price
100	\$59	100	\$88

To order standard packs, add a -M after the item #.



Case using a technique known as Deep Margin Elevation (DME) or cervical margin relocation that has been evolving since at least the 1990s

Dr. Robert Rosenfeld (Westwood, New Jersey) case utilizing the Garrison® Reel Matrix™ margin elevation band.



Figure 1. Removal of the amalgam and excavation revealed that the decay extended quite far subgingivally on both the mesial and distal.



Figure 2. After placing a self etching primer and adhesive, I applied several small increments of a bulk fill flowable material, enough to rise above the gingiva.




Figure 3. I then did a freehand layered buildup of the remainder of the restoration, with the exception of the mesial box, using a light-cured core buildup material.




Figures 4 and 5. At this point, I removed the Margin Elevation band and placed a sectional matrix band, wedge and separating ring on the mesial (Composi-Tight® 3D Fusion™, Garrison® Dental) and completed the buildup to create good contour and contact until such time that the tooth is prepared for a crown.





Figure 6. After trimming the slight excess material on the distal with a finishing carbide and adjusting the occlusion, here is the finished restoration.

How PerForm™ works:



STEP 1:
Place the first composite increment. Press PerForm™ into the uncured composite increment and either push (distal; fig. 1A) or pull (mesial; fig. 1B).



STEP 2: Place the curing light directly on the Fresnel lens and cure.

PerForm

Proximal Contact Instruments

PerForm puts you in control.

- Improves contacts with any matrix system
- Control the height and breadth of contacts
- Increase the tightness of contacts
- Unique lens controls curing light diffusion
- Steam autoclave only

PerForm™

Item	Kit	Qty.	Price
CFK01	PerForm (1/ea. size)	2/pk	\$126
CF100	1.1 mm Small Instruments	2/pk	\$126
CF200	2 mm Large Instruments	2/pk	\$126



"To ensure a tight contact a bulk fill flowable was placed in one proximal box at a time and the mylar held against the adjacent tooth with a PerForm™ proximal contact instrument.

This ensures that the contact is tight and well-formed while light curing through it."

*Dr. Jack D. Griffin, Jr., DMD, MAGD
Weldon Spring, MO*



Photo courtesy of David S. Hornbrook, DDS, FAACD International Lecturer and Private Practice, San Diego, CA


Curved Form/tapered: Provides the optimum curvature and band height for virtually every anterior restoration.

Blue Tint: Provides contrast between matrix and tooth structure without compromising composite polymerization.

Our Most Popular Anterior Matrix

BlueView™ VariStrip™

Contoured Anterior Matrices



The Blue View™ VariStrip™ is a revolutionary contoured anterior matrix providing the optimum curvature and band height for virtually every anterior restoration. The 0.0020" thin plastic anatomical strip is tapered from one end to the other. You can place the strip interproximally and then slide it until the tooth height is perfectly matched. The pre-contouring easily recreates the occluso-gingival anatomy, avoiding flat embrasures.

Blue View™ VariStrip™

Item	Kit	Qty.	Price
VS01	VariStrip	100	\$117

BlueView™

Cervical Matrices

"I love the Blue View Cervical Matrix! It works so well. My assistants love watching me use this system. The patients love the excitement in our faces. It makes the process easy and fast. Gotta love that! We are also hooked on using the ReelMatrix. Thanks so much for making great products."

Cynthia Brattesani, DDS,
San Francisco, CA, USA



Flexible

- adapts to individual anatomy for improved contours
- use it to directly compact the composite into the restoration
- thin edges allow excess composite to flow out for improved anatomy

Blue Tint


- provides contrast between matrix and tooth structure
- does not compromise composite polymerization

Convenient!

- protects the restoration from contamination
- eliminates time spent hand sculpting
- eliminates air-inhibited layer
- reduces finishing time
- 7 sizes for a perfect fit
- Gingival Retractor included!

KCMK1 INCLUDES







- 275 Assorted Matrices
- 1 Matrix placement Instrument
- 1 Gingival Retractor



MATRICES INCLUDED:

- Anterior/Premolar
KCMS60 - 50 pcs
- Molar
KCMS70 - 50 pcs
- KCMS10 - 35 pcs
- KCMS20 - 35 pcs
- KCMS30 - 35 pcs
- KCMS40 - 35 pcs
- KCMS50 - 35 pcs

Blue View™ Cervical Matrices

Item	Description	Price
KCMK1	Assortment Kit w/Instruments	\$188
Refills		
KCMS10	 Refill 100 pack	\$104
KCMS20	 Refill 100 pack	\$104
KCMS30	 Refill 100 pack	\$104
KCMS40	 Refill 100 pack	\$104
KCMS60	 Refill 100 pack	\$104
KCMS70	 Refill 100 pack	\$104
KCMINST1	Placement Instrument	\$21

Composi-Tight® Matrix Forceps

Precise Placement, Unbeatable Grip.



Check out the GRIP! Effortlessly remove the matrix bands that otherwise would be a struggle.



Precise Placement: The tips have a 75-degree angle with an occlusal-gingival orientation, providing a more appropriate placement angle than cotton forceps. The petite, carbide tips allow the matrix to be grasped at its upper edge to prevent distortion of the interproximal area of the band.



Unbeatable Grip: The tungsten carbide tips provide exceptional grip and tactile control for tough access cases.

Matrix Forceps

BFA-1 Angled tip for occlusal placement \$160



Precision Instrumentation completes the NiTin™ system!

Low profile tips simplify stacking on MODs with NiTin's smaller rings. Recommended for NiTin™ rings only.



NiTIn™ Precision Instruments

NTG-MRDF-100	Matrix Forceps	\$79
NTG-DF6-8	Surgical Tissue Forceps	\$27

This surgical tissue forceps with its deep interdigitating teeth provide ample grip and tactile control for both band and wedge placement.

toll free 888.437.0032 | www.garrisondental.com

Composi-Tight® 3DFusion™

Composi-Tight® 3D Fusion™ includes dentistry's premier heavy-duty stainless steel ring placement forceps. Regardless of your hand size or strength, Composi-Tight® 3D Fusion™ forceps will give you incredible control for the precise placement of super-strong separator rings over even the widest molars.


Composi-Tight® 3D Fusion™ Ring Placement Forceps
FXP01 \$102

Composi-Tight® 3D Fusion™ forceps will work with virtually any ring from any manufacturer.



Ring Placement Forceps

Holds all your sectional matrix system rings.



Ring Placement Forceps
AUMRDF-100 \$85

Universal Forceps
for placing rubber dam clamps and rings
UNVF-100 \$85



AUMRDF-100



UNVF-100

Composite Instruments Nonstick titanium-nitride coated

Our multi-function instruments have double condenser configurations and carving blades on opposite ends. Now you can place, pack, shape and carve your way to perfectly contoured composite restorations without ever having to switch instruments!

- Large-diameter, lightweight handles for improved ergonomics.
- Provides non-stick surface for placement and shaping of composites.
- 40% smoother than stainless steel for a non-porous, non-stick surface.
- Withstands diverse cleaning/sterilization methods, including ultrasonic cleaning.



"Skippy Special"

TN008

Item	Description	Price
TN008	"Skippy Special" 2 condensers, ball burnisher, blade, acorn	\$108



5-in-1

TN009

Item	Description	Price
TN009	5-in-1 Multi-Function Instrument 2 blades, 2 condensers, acorn	\$108

Item	Description	Price
TN001	Medium Blades Placement of composite material Gingival retraction to sculpt under free margins Burnishing matrix bands	\$46



TN002	Medium Blade with Medium Plugger Contouring small Class I, II, V Packing composite material into conservative preparations Burnishing matrix bands	\$46
-------	--	------



TN003	Medium Blade with Large Plugger Contouring large Class I, II, V Packing composite material into large preparations Convex plugger tip minimizes pullback of composite material	\$46
-------	--	------



Item	Description	Price
TN004	Extra Thin Blades Reach under margins to hone and sculpt Sculpting interproximally Packing cord in gingival crevice	\$49



TN005	Long Extra Thin Blades Allows additional access for hard-to-reach posterior contouring and shaping Sculpting interproximally Packing cord in gingival crevice	\$49
-------	---	------











TN006	Small/Medium Multi-Function Plugger Placement of composite material in smaller Class II and III Packing material toward gingival seats Carving occlusal anatomy with acorn carver	\$75
-------	---	------



Matrix Selection Guide

Being the Matrix Company, one of the most common questions we get is “which matrix do I use where?” While personal preference will drive much of the decision process, below is a chart that can provide a good starting point.

Procedure	Recommended Materials	Page	Special Notes
Distal of canine, short or malpositioned 	Composi-Tight® 3D Fusion™ Sectional Matrix Kit	4–5	MOST POPULAR! Separator rings designed to hold firmly on tough cases.
	Reel Matrix™—Metal or Blue View	30–31	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Slick Bands™ Tofflemire-Style Matrices	32–33	Traditional Tofflemire-style with non-stick coating.
Conservative Class II 	Composi-Tight® 3D Fusion™ Sectional Matrix Kit	4–5	MOST POPULAR! Superior contacts with minimal flash.
	Reel Matrix™—Metal or Blue View	30–31	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Slick Bands™ Tofflemire-Style Matrices	32–33	Traditional Tofflemire-style with non-stick coating.
Slightly wide and/or deep 	Composi-Tight® 3D Fusion™ Sectional Matrix Kit	4–5	MOST POPULAR! Superior contacts with minimal flash.
	Reel Matrix™—Metal or Blue View	30–31	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Slick Bands™ Tofflemire-Style Matrices	32–33	Traditional Tofflemire-style with non-stick coating.

Procedure	Recommended Materials	Page	Special Notes
Large Class II with missing cusps 	Composi-Tight® 3D Fusion™ Sectional Matrix Kit—Green Wide Prep Ring	4–5	EXCLUSIVE! Dentistry's first Wide Prep ring greatly simplifies these more challenging restorations.
	Reel Matrix™—Metal or Blue View	30–31	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Slick Bands™ Tofflemire-Style Matrices	32–33	Traditional Tofflemire-style with non-stick coating.
Small to medium size Class III 	Fusion™ Anterior Matrix System	7	Bands and wedges with ideal gingival-incisal curvature for many types of anterior lesions.
	VariStrip™ Contoured Anterior Matrices	36	Pre-contoured in three dimensions AND tapered to provide the ideal band height and contour.
Large Class III or Class IV 	Fusion™ Anterior Matrix System	7	Bands and wedges with ideal gingival-incisal curvature for many types of anterior lesions.
	VariStrip™ Contoured Anterior Matrices	36	Pre-contoured in three dimensions AND tapered to provide the ideal band height and contour.
Class V 	Blue View™ Cervical Matrix	37	Contours and directly compresses resin while eliminating air-inhibited layer.
Crown build-ups, post and core 	Reel Matrix™	30–31	Best choice. Margin Elevation band available for very deep buildups. Retainerless convenience.
	Pinch™ Matrices—Metal or Blue View	34	Very fast; may not provide as tight a cervical seal as other matrices.
	Slick Bands™ Tofflemire-Style Non-Stick Matrices	32–33	Margin Elevation band available for very deep buildups if you prefer a traditional matrix retainer.

Order Form

Fax: 616.842.2430

Item Number	Description	Quantity	Price/Unit	Total Price
		Shipping		
		Sales Tax		
		Total Price		

Doctor's Name _____

Contact Name _____

Address _____

City _____ State _____ Zip _____

Phone _____ Fax _____

Email _____

Payment Method

Credit Card: Visa MasterCard Am. Express Discover C.O.D. (add \$13.00)

Credit Card # _____ CC Ref # _____

Exp. Date _____ Signature _____

Credit Card Billing Information

Address _____

City _____ State _____ Zip _____

Call for additional express shipping rates and delivery. All prices subject to change without notice. Appropriate shipping charges and sales tax will be added. Special offers not valid with any other advertised specials.

toll free 888.437.0032 | garrisondental.com

Customer Service Information

Call **888.437.0032 Toll Free**
(USA & Canada)

616.842.2244

Fax 616.842.2430

Email gds@garrisondental.com

Online **garrisondental.com**
Look for our Learning Center with clinical tips, case studies and technique guides.

Mail 150 DeWitt Lane
Spring Lake, MI 49456
USA

Shipping Rates

Regular UPS	2nd Day Air	Next Day Air	C.O.D.
\$9.75	\$19.00	\$34.00	+\$13.00

Express orders must be placed by 4pm EST to ship same day.

Shipping charges subject to change without notice.

Payment Methods

VISA, MasterCard, American Express,
Discover, Check Prepay or C.O.D.

Prices are subject to change without notice.

Distribution

Garrison is a global leader for practical and innovative dental products and accessories, with distributors worldwide.

We are pleased to announce that your favorite Garrison products are available through several full-service dealers in the USA.

To find a dealer, contact our headquarters in Spring Lake, Michigan. Phone: 616.842.2244 or Email: gds@garrisondental.com



Order Toll-Free • Order by Fax • Order Online • Order by Mail

Our Guarantee – It's Simple!

Satisfaction = 100%

Everything you order from Garrison® Dental Solutions is backed by our risk-free six-month money-back **Satisfaction Guarantee**. If you are not satisfied with your product, return the unused portion and we will issue a refund.

Quality = 100%

All Garrison products are covered by our **Quality Guarantee**. We will repair or replace any Garrison® product that is shown to be defective in materials or workmanship based on the individual product's anticipated lifetime. Please visit our website for details on our various products, garrisondental.com.

Garrison[®]
Dental Solutions

150 DeWitt Lane,
Spring Lake, MI 49456 USA
toll free 888.437.0032
www.garrisondental.com

RLCAT122